


Tripod Stand

EP-3572



Equipment

ITEM	ID	ITEM	ID
Vertical Rail	1	Center Leg	3
Left Leg	2	Right Leg	4


Introduction

The EP-3572 Tripod Stand for Essential Physics 3rd Ed consists of a vertical rail and three legs. The vertical rail has “T-slots” on four sides and there are a total of five nut sliders already installed in two of the “T-slots”. The nut divider is designed to hold 1/4-20 threaded rods and bolts such as the ones on other components in the Essential Physics Trays, Kits and Modules.

Assembly

Each leg has a knob and a threaded bolt at the tall end. The head of the threaded bolt slides into the “T-slot” on the side of the vertical rail.

Turn the knob counter-clockwise to loosen the threaded bolt so that the bolt head is two or three millimeters away from the leg. Align the head of the threaded bolt with a “T-slot” on the bottom end of the vertical rail. Slide the head of the bolt into the “T-slot”. Continue to slide the leg onto the rail until the lip at the bottom of the leg touches the bottom of the rail.



Turn the knob clockwise to tighten the leg onto the vertical rail. Repeat the process for the other legs.

Operation

Nut Sliders

The nut sliders are in the “T-slots” of the vertical rail and are held in place by friction. They slide up or down along the vertical rail.

Connect a compatible threaded rod to a nut slider. Tighten the rod or bolt to hold the nut slider in place. Loosen the rod or bolt to move the component and nut slider up or down the rail.

Remove one or more legs in order to move nut sliders from one “T-slot” to another. Slide the nut slider to the bottom of the vertical rail to remove it. Carefully align it with the bottom of the other “T-slot” and slide it up the vertical rail.

Technical Support

For assistance with any PASCO product, contact PASCO at:

Address: PASCO scientific
10101 Foothills Blvd.
Roseville, CA 95747-7100

Phone: 916-462-8384 (worldwide)
800-772-8700 (U.S)

Email: techsupp@pasco.com

Web: www.pasco.com/support

Limited Warranty

For a description of the product warranty, see the PASCO catalog.

Copyright

The PASCO scientific Instruction Sheet is copyrighted with all rights reserved. Permission is granted to non-profit educational institutions for reproduction of any part of this manual, providing the reproductions are used only in their laboratories and classrooms, and are not sold for profit. Reproduction under any other circumstances, without the written consent of PASCO scientific, is prohibited. Rev: 01/18

Trademarks

PASCO, PASCO scientific, PASCO Capstone and SPARKvue are trademarks or registered trademarks of PASCO scientific, in the United States and/or in other countries. All other brands, products, or service names are or may be trademarks or service marks of, and are used to identify, products or services of, their respective owners. For more information visit www.pasco.com/legal.